Stock Code: 688068.SH

Novel Coronavirus 2019-nCoV Antigen Test (Colloidal Gold) Summary Data

Beijing Hotgen Biotech Co., Ltd.

Explanation of The Export License

To whom it may concern

According to the No.12 Government Notice published by MOC, GAC and NMPA, CCCMHPIE announced the list of companies which were permitted to export novel coronavirus diagnostic products, Hotgen is in the list as aforesaid.

The screenshot shows the website of the China Pharmaceutical and Health Products Import and Export Association (CCCMHPIE). The page title is "中国医药保健品进出口商会" (China Pharmaceutical and Health Products Import and Export Association). The search bar contains the text "91110115777090586H". The table below lists the companies permitted to export novel coronavirus diagnostic products.

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
Novel Coronavirus 2019-nCoV Antigen Test (Colloidal Gold)

Product Features

- High Accuracy, Specificity and Sensitivity
- No need instrument, get results in 15 minutes
- Room temperature storage
- Sample : Human Anterior Nares Swab
- Detect the presence of viral proteins
- Identify acute or early infection

Clinical Performance

(Disease Course 5-7 Days)

Sensitivity: 96.30%; Specificity: 99.13%; Accuracy: 97.76%.

Novel Coronavirus 2019-nCoV Antigen Test(Colloidal Gold)

Specimen Requirements

1 Sample collection

Gently insert the entire soft tip of the swab into one nostril for 1.5cm until you feel a bit of resistance.


Using medium pressure, rub the swab slowly in a circular motion around the inside wall of your nostril 4 times for a total time of 15 seconds.

Repeat the same process with the same swab in the other nostril.

2 Sample treatment

The swab after sampling is soaked below the liquid level of the sample extraction buffer. Rotate and press 3 times. The swab soaking time is not less than 15s.

The swab head is pressed, then take out the swab and tighten the sampling tube.

3 Sample preservation

The treated sample should be tested within 1h.

Test Procedure

Place the test cassette, sample extraction buffer at room temperature for 15-30 minutes, and equilibrate to room temperature (10-30 °C).

Open the aluminum foil pouch of the test cassette, place the test cassette on a flat surface.




Add 4 drops of the treated sample into the sample well of the test cassette. (In case of chromatographic abnormalities, extra add 1-2 drops of the treated sample accordingly.) Incubate at 10-30 °C for 15 minutes.

Observe the results after incubate at 10-30 °C for 15 minutes. The result after 30 minutes is invalid.

Interpretation of result

Clinical Performance

This study enrolled a total of 223 human anterior nasal swab specimens, of which 108 were positives for RT-PCR Test, 115 negatives for RT-PCR Test; As for data collection of the corresponding RT-PCR Test results.

Assessment system Novel Coronavirus 2019-nCoV Antigen Test (Colloidal Gold)	Reference system (clinical diagnostic results)		
	Positive (+)	Negative (-)	Total
Positive(+)	104	1	105
Negative(-)	4	114	118
Total	108	115	223

Sensitivity: 96.30%; Specificity: 99.13 %; Accuracy: 97.76%.

Product information

Product name	Test samples	Specifications	Storage conditions
Novel Coronavirus 2019-nCoV Antigen Test (Colloidal Gold)	Human Anterior Nares Swab	1T/kit, 5T/kit, 20T/kit, 40T/kit	4-30 °C

Beijing Hotgen Biotech Co.,Ltd

Add.: 9 building, No.9 Tianfu Street, Daxing District, Beijing, 102600, P.R. China

Website: www.hotgen.com.cn

Declaration of Conformity

Manufacturer:

Name: Beijing Hotgen Biotech Co.,Ltd

Address: 9th building, No.9 Tianfu Street, Biomedical Base,Daxing District, Beijing,
102600, P.R.China

European Representative:

MedNet GmbH

Borkstrasse 10,48163 Muenster,Germany

Product Name:

Novel Coronavirus 2019-nCoV Antigen Test (Colloidal Gold)

Novel Coronavirus 2019-nCoV Antigen Test (Up-converting Phosphor Immunochromatographic
Technology)

Classification : **Others of ANNEX II of IVDD**

Conformity Assessment Route: **Annex III**

We herewith declare that the above mentioned products meet the transposition into national law, the provisions of the following EC Council Directives and Standards. All supporting documentations are retained under the premises of the manufacturer.

General applicable directives:

In Vitro Diagnostic Medical Devices DIRECTIVE 98/79/EC

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Formblatt für In-vitro-Diagnostika / Form for In Vitro Diagnostic Medical Devices

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Instructions for Use

FREQUENTLY ASKED QUESTIONS

- When can I test myself?
- You can always test yourself whether you have symptoms or not. Please note that the test results are a snapshot that is valid for this point in time. Tests should therefore be repeated according to the regulations of the responsible authorities.
- What should I pay attention to in order to obtain the most exact test result possible?
- Always follow the instructions for use exactly. Perform the test immediately after collecting the sample. Dispense the drops from the test tube only into the designated well of the test cassette. Dispense four drops from the sample tube. Too many or too few drops can lead to an incorrect or invalid test result.
- The test strip is very discolored. What is the reason or what am I doing wrong?

The reason for a clearly visible discoloration of the test strip is that too large a quantity of drops has been dispensed from the sample tube into the test cassette well. The indicator strip can only hold a limited amount of liquid. If the control line does not appear or the test strip is very discolored, please repeat the test with a new test kit according to the instructions for use.

- What should I do if I took the test but didn't see a control line?
- In this case, the test results to be considered invalid. Please repeat the test with a new test kit according to the instructions for use.
- I am unsure of the interpretation of the results. What should I do?
- If you cannot clearly determine the result of the test, contact the nearest medical facility applying the regulations of your local authority.
- My results are positive. What should I do?
- A horizontal colored line is visible in the control area (C) as well as in the test area (T), your results are positive and you should immediately contact the medical facility in accordance with the requirements of your local authorities. Your test result may be checked and the next steps will be explained to you.
- My results are negative. What should I do?
- If only a horizontal colored line is visible in the control area (C), this may mean that you are negative or that the viral load is too low to be recognized by the test. If you experience symptoms such as headaches, migraines, fever, loss of sense of smell and taste, contact the nearest medical facility applying the regulations of your local authority. In addition, you can repeat the test with a new test kit.
- Can this test cassette be reused or used by multiple people?
- This test cassette is for one-time use and cannot be reused or used by multiple people.

MODEL NUMBER

Model A

SPECIFICATIONS

1T/kit, 5T/kit, 20T/kit, 25T/kit, 40T/kit, 50T/kit.

INTENDED USE

This kit is used for in vitro qualitative determination of SARS-CoV-2 antigens in human anterior nasal swab samples. It can be used for rapid investigation of suspected COVID-19 cases, and can also be used as a reconfirmation method for nucleic acid detection in discharged cases. A positive test result indicates that the sample contains SARS-CoV-2 antigen. A negative test result does not rule out the possibility of infection. This kit is for home use by laymen in a non-laboratory setting (such as persons' home or certain non-traditional sites such as offices, sporting events, airports, schools etc.). The test results of this kit are for clinical reference only. It is recommended to conduct a comprehensive analysis of the condition based on

the patient's clinical manifestations and other laboratory tests.

COMPONENTS

1. SARS-CoV-2 Antigen Test Cassette
2. Sample extraction buffer
3. Disposable virus sampling swab
4. Biohazard specimen bag

Note: Components of different batches cannot be mixed.

SPECIMEN REQUIREMENTS

1. Sample collection

- Gently insert the entire soft tip of the swab into one nostril for 1.5cm until you feel a bit of resistance.
- Using medium pressure, rub the swab slowly in a circular motion around the inside wall of your nostril 4 times for a total time of 15 seconds.
- Repeat the same process with the same swab in the other nostril.

2. Sample treatment

- The swab after sampling is soaked below the liquid level of the sample extraction buffer. Rotate and press 3 times. The swab soaking time is not less than 15 seconds.
- The swab head is pressed, then take out the swab and tighten the sampling tube.

3. Sample preservation: The treated sample should be tested within 1h.

TEST PROCEDURE

1. Place the test cassette, sample extraction buffer at room temperature for 15~30 minutes, and equilibrate to room temperature (10~30°C).

2. Open the aluminum foil pouch of the test cassette, place the test cassette on a flat surface.

3. Add 4 drops of the treated sample into the sample well of the test cassette. (In case of chromatographic abnormalities, add an extra of 1~2 drops of the treated sample accordingly). Incubate at 10~30°C for 15 minutes.

4. Observe the results after incubating at 10~30°C for 15 minutes. The result

obtained after 30 minutes is invalid.

15 min.

DISPOSAL THE SAMPLE AND CLEAN-UP

- Place the test cassette, sample extraction buffer and disposable virus sampling swab in the biohazard specimen bag and seal the bag.
- Throw away the remaining sample kit items.
- Re-apply hand sanitizer.

INTERPRETATION OF RESULT

Positive: Two color bands appear in the observation window, that is, a red or magenta line appears at the position of the quality control line (C line) and the detection line (T line) (as shown in result 1), indicating the test result of SARS-CoV-2 antigens in the sample is positive.

Negative: A red or magenta line appears at the position of the quality control line (C line) in the observation window, and no line appears at the position of the test line (T line) (as shown in the result 2), indicating the test result of SARS-CoV-2 antigens in the sample is negative or the concentration is below the limit of detection of the kit.

Invalid: No line appears in the position of the quality control line (line C) in the observation window (as shown in result 3), indicating that the test is invalid, and the sample should be recollect and re-tested.

Result 1: Positive **Result 2: Negative** **Result 3: Invalid**

PRINCIPLE OF THE ASSAY

This kit is based on the colloidal gold immunochromatographic technology, and uses the double antibody sandwich method to detect N protein of SARS-CoV-2 antigen in human anterior nasal swab samples. The detection line (T line) of the SARS-CoV-2 antigen test cassette was coated with SARS-CoV-2 antibody, and the quality control line (C line) was coated with sheep anti-mouse antibody. During the test, the sample is dropped into the test cassette and the liquid is chromatographed upward under the capillary effect. The SARS-CoV-2 antigen in the sample first binds to the colloidal gold-labelled SARS-CoV-2 antibody to form a solid phase SARS-CoV-2 antibody-SARS-CoV-2 antigen-labelled SARS-CoV-2 antibody-colloidal gold complex at the T line position, and form a solid phase sheep anti-mouse-labelled SARS-CoV-2 antibody-colloidal gold complex was formed at the C line position. After the test is completed, observe the colloidal

gold color reaction of T line and C line to determine results of SARS-CoV-2 antigen in human anterior nasal swab samples.

STORAGE AND SHELF LIFE

- The kit should be stored at 4~ 30°C, the shelflife is setfor 18 months.
- After the foil bag is opened, it should be used within 30 minutes (temperature 10~30°C, humidity ≤70%).
- The sample extraction buffer should be used within 18 months after opening (temperature 10~30°C, humidity ≤70%).

See label for manufacture date and expiration date.

LIMITATIONS

- The test result of this kit is not the only confirmation indicator of clinical indications. The infection should be confirmed by a specialist along with other laboratory results, clinical symptoms, epidemiology, and additional clinical data.
- In the early stages of infection, low levels of antigen expression can result in negative results.
- The sample test results are related to the quality of sample collection, processing, transportation and storage. Any errors may lead to inaccurate results. If cross-contamination is not controlled during sample processing, false positive results may occur.

PERFORMANCE CHARACTERISTICS

- Limit of Detection (LoD)
Novel Coronavirus 2019-nCoV Antigen Test (Colloidal Gold) has been confirmed can detect SARS-CoV-2 at $2.5 \times 10^{2.2}$ TCID₅₀/mL, which was collected from a confirmed COVID-19 patient in China.
- Study on Exogenous/Endogenous Interference Substances:
The potential interfering substances listed below do not interfere.

(1) Exogenous factor

No.	Exogenous factor	Interfering substances	Test conc.
1	Nasal sprays or drops	Phenylephrine Oxymetazoline	128 µg/mL 128 µg/mL
2		Saline Nasal Spray 10%	10% (v/v)
3		Dexamethasone	2 µg/mL
4	Nasal corticosteroids	Flunisolide	0.2 µg/mL
5		Triamcinolone acetate	0.2 µg/mL
6		Mometasone	0.5 µg/mL
7		Stripsils (flurbiprofen 8.75mg)	5% (w/v, 50mg/mL)
8	Throat lozenges	Throat candy	5% (w/v, 50mg/mL)
9		Anbesol (Benzocaine 20%)	5% (v/v)
10	Oral anaesthetic	α-Interferon-2b	0.01 µg/mL
11		Zanamivir (Influenza)	2 µg/mL
12		Ribavirin (HCV)	0.2 µg/mL
13		Osetamivir (Influenza)	2 µg/mL
14	Anti-viral drugs	Peramivir (Influenza)	60 µg/mL
15		Loximavir (HIV)	80 µg/mL
16		Ritonavir (HIV)	20 µg/mL
17		Atazanavir (HIV)	40 µg/mL
18		Levofloxacin Tablets	40 µg/mL
19		Azithromycin	200 µg/mL
20	Antibiotic	Ceftriaxone	800 µg/mL
21		Metoprolol	100 µg/mL
22	Antibacterial, systemic	Tobramycin	128 µg/mL
23		Mucin: bovine submaxillary gland, type	100 µg/mL
24	Other	Biotin	100 µg/mL
25			

(2) Endogenous factor

No.	Endogenous factor	Interfering substances	Test conc.
1	Autoimmune disease	Human anti-mouse antibody, HAMA	800 ng/mL

2	Serum protein	Whole blood (human), EDTA anticoagulated	10% (w/w)
---	---------------	--	-----------

3. Cross-Reactivity & Microbial interference:

There is no cross-reaction and no interference with the potentially cross-reactive microorganisms listed below.

No.	Crossing reacting substance	Strain	Concentration of cross reacting substance
1	Human Coronavirus	HKU1	2×10^5 TCID ₅₀ /mL
2		229E	2×10^5 TCID ₅₀ /mL
3		OC43	2×10^5 TCID ₅₀ /mL
4		NL63	2×10^5 TCID ₅₀ /mL
5		SARS	2×10^5 TCID ₅₀ /mL
6		MERS	2×10^5 TCID ₅₀ /mL
7		Type 1	2×10^5 TCID ₅₀ /mL
8		Type 2	2×10^5 TCID ₅₀ /mL
9		Type 3	2×10^5 TCID ₅₀ /mL
10	Adenovirus	Type 4	2×10^5 TCID ₅₀ /mL
11		Type 5	2×10^5 TCID ₅₀ /mL
12		Type 7	2×10^5 TCID ₅₀ /mL
13		Type 55	2×10^5 TCID ₅₀ /mL
14	Human Metapneumovirus (hMPV)	hMPV 3, Type 81 / Perth-2002	2×10^5 TCID ₅₀ /mL
15		hMPV 16 Type A1 / IA1D-2003	2×10^5 TCID ₅₀ /mL
16		Type 1	2×10^5 TCID ₅₀ /mL
17		Type 2	2×10^5 TCID ₅₀ /mL
18	Parainfluenza virus	Type 3	2×10^5 TCID ₅₀ /mL
19		Type 4A	2×10^5 TCID ₅₀ /mL
20		H1N1	2×10^5 TCID ₅₀ /mL
21	Influenza A	H2N2	2×10^5 TCID ₅₀ /mL
22		H5N1	2×10^5 TCID ₅₀ /mL
23		H7N9	2×10^5 TCID ₅₀ /mL
24		Yamagata	2×10^5 TCID ₅₀ /mL
25		Victoria	2×10^5 TCID ₅₀ /mL
26		Type 6B	2×10^5 TCID ₅₀ /mL
27	Enterovirus	09/2014 isolate 4	2×10^5 TCID ₅₀ /mL
28	Respiratory syncytial virus	Type A	2×10^5 TCID ₅₀ /mL
29		A16	2×10^5 TCID ₅₀ /mL
30		Type B42	2×10^5 TCID ₅₀ /mL
31	Rhinovirus	TWAR strain TW-183	5×10^6 CFU/mL
32	Chlamydia pneumoniae	NCTC 4560	5×10^6 CFU/mL
33	Haemophilus influenzae	Bloomington-2	5×10^6 CFU/mL
34		Los Angeles-1	5×10^6 CFU/mL
35	Legionella pneumophila	82A3105	5×10^6 CFU/mL
36		K	5×10^6 CFU/mL
37		Erdman	5×10^6 CFU/mL
38	Mycobacterium tuberculosis	HN878	5×10^6 CFU/mL
39		CDC1551	5×10^6 CFU/mL
40		H37Rv	5×10^6 CFU/mL
41		4752-98 (DJ)68-17	5×10^6 CFU/mL
42	Streptococcus pneumoniae	178 [Poland 23F-16]	5×10^6 CFU/mL
43		262 [CIP 104340]	5×10^6 CFU/mL
44		Slovakia 14-10 [29055]	5×10^6 CFU/mL
45		Typing strain T1 [NCIB 11843, 5F 130]	5×10^6 CFU/mL
46	Streptococcus pyogenes	NCP 13671	5×10^6 CFU/mL
47	Bordetella pertussis	Mutant 22	5×10^6 CFU/mL
48	Mycoplasma		

49	Pneumoniae	PH strain of Eaton Agent [NCTC 10119]	5×10^6 CFU/mL
50		M129-B7	5×10^6 CFU/mL
51	Pneumocystis jirovecii (PJP)	N/A	N/A
52	Pooled human nasal wash	N/A	N/A
53	Candida albicans	3147	5×10^6 CFU/mL
54	Pseudomonas aeruginosa	R. Hugh 813	5×10^6 CFU/mL
55	Staphylococcus epidermidis	FDA strain PCI 1200	5×10^6 CFU/mL
56	Streptococcus salivarius	S21B [FO 13956]	5×10^6 CFU/mL

4. Hook Effect:
There is no hook effect at $1.0 \times 10^{2.2}$ TCID₅₀/mL of SARS-CoV-2 isolated from a SARS-CoV-2 confirmed patient in China.

5. Clinical Performance:

Clinical performance of Novel Coronavirus 2019-nCoV Antigen Test (Colloidal Gold) has been determined by testing 108 positive and 115 negative specimens for SARS-CoV-2 antigen (Ag). The sensitivity is 96.30% (95% CI: 90.79-98.98%), and the specificity is 99.13% (95% CI: 95.25-99.98%).

Novel Coronavirus 2019-nCoV Antigen Test (Colloidal Gold) Res ults	PCR Test Results		Total
	Positive	Negative	
Positive	104	1	105
Negative	4	114	118
Total	108	115	223

Sensitivity	Specificity	Overall Percentage Agreement
96.30% (90.79%-98.98%)	99.13% (95.25%-99.98%)	97.76% (94.85%-99.27%)

PRECAUTIONS

- This kit is for in vitro diagnostic use only. Please read this instruction carefully before the test.
- Please use the swab and sample extraction buffer provided in this kit, and do not replace the sample extract in this kit with components in other kits.
- Operations should strictly follow the instructions.
- Positive and negative predictive values are highly dependent on the prevalence. When the prevalence of the disease is low and SARS-CoV-2 has little/no activity, a positive test result is more likely to represent a false positive result; when the prevalence of the disease is high, false negative test results are more likely.
- Compared with a RT-PCR SARS-CoV-2 assay, this test is less sensitive when used to detect patient samples within the first five days of the onset of symptoms.
- The test cassette must be used within 30 minutes after opening (temperature 10~30°C, humidity ≤70%), it should be used immediately after opening at 30°C, and the unused test cassette must be sealed and dryly stored.
- Waste or excess samples produced during testing should be inactivated according to regulations on infectious agents.

EXPLANATION FOR IDENTIFICATION

	Use by date		Batch		Consult instruction for use
	Content Sufficient For n Tests		Temperature limitation		Catalog Number
	Manufacturing date		Caution		Do not reuse

	CE Marking – IVDD 98/79/EC		Authorized representative in the European Community		Manufacturer
	For In Vitro Diagnostic Use		Keep away from sunlight		Keep dry
	For self-testing	/	/	/	/

Beijing Hotgen Biotech Co., Ltd.
 9th Building, No. 9 Tianfu Street, Biomedical Base,
 Daxing District, Beijing, 102600, P.R. China.

MedNet GmbH
 Borkstrasse 10, 48163 Muenster, Germany

APPROVAL DATE AND REVISION DATE OF THE INSTRUCTION

Approved on February, 2021;

Version number: V. 2021-02.01 [Eng.]